

The Power of Light!

CR-2

Digital Retinal Camera



- Non-Mydriatic Color Imaging
- Compact & Lightweight (33 lbs)
- Minimum Pupil Size at 3.3 mm (Small Pupil Mode)
- Low Flash Intensity using White LED Flash
- 45° Angle with 2X Digital Magnification (30-degree)
- Equipped with 15.1 MP Canon EOS Digital Camera for CR-2
- Auto Exposure Mode Adjusts Brightness of IR View and Acquired Image
- Retinal Imaging Control Software (RICS) Provides DICOM and JPEG Output
- RICS Provides Digital Output of Cobalt, Red-Free, RGB and Cup/Disc Ratio, Stereo Pair
- Anterior Segment Lens Compensator Supports External Eye Photography

Canon