


**Now Integrates with the Huvitz Digital Refraction System (see page 3)**

# Sirius™ Visual Acuity System


The Sirius Visual Acuity System offers a complete set of visual acuity tests in one easy-to-use product. Choose from a wide range of acuity tests that can be randomized to prevent memorization and improve acuity testing accuracy.

All optotypes may displayed as:

- Single character mode (crowding bars on/off)
- Single line mode (horizontal and vertical)
- Multi line mode
- Single column mode (speed acuity testing)



Red Green Mask



Tumbling E's



Contrast Sensitivity



Allen Preschool



## Acuity Charts

- Optotypes
- Worth 4 dot
- Contrast sensitivity
- Astigmatism
- Circular astigmatism
- Concentric circles
- Red green balance
- Fixation points

## Optotypes

Measurement in Snellen and ETDRS

- HOTV
- Allen Preschool
- Tumbling E
- Landolt C
- Numbers



A default video and image gallery is included as a standard feature. This can be used as a child's fixation device and for educational purposes. In addition, the galleries can be customized to your specific needs to include your own unique images and video presentations.

Infrared remote control for convenience and ease of operation.