


## ∴ HEINE OMEGA 600 Breath Shield

Designed to cover the mouth-nose-chin area. For use with the HEINE OMEGA 600 only.


Easy to attach.



1. Align the breath shield with OMEGA 600 optics.



2. Press Breath Shield onto OMEGA 600 until securely attached.



3. Once firmly attached, the shield should not move easily.



4. Take the Breath Shield off for cleaning and disinfection.

More information: [heine.com/en/cleaning-breath-shield-omega-600](https://heine.com/en/cleaning-breath-shield-omega-600)

### Article numbers

- ∴ HEINE OMEGA 600 Breath Shield – 1 x: **C-000.33.018**
- ∴ HEINE OMEGA 600 Breath Shield – 3 x: **C-000.33.019**