

# Advanced Features, Creating Better Value

## Stimuli™ Advanced Version 4.0



Stimuli™ is a refraction device that replaces the old technology, the automated projector. Stimuli™ eliminates the "memore effect" by bringing virtually unlimited charts to the patient in a random manner. All functions are at the practitioners' fingertips.



Since its introduction in 2000, the Stimuli™ Visual Acuity System has been a very reliable tool in the hands of physicians in clinics as well as private practice. The new 21.5" iMac slim line makes Stimuli™ Advanced an even more desirable and versatile diagnostic tool.

### **Stimuli™ Advanced Version 4.0 Includes:**

- > 21.5" Apple iMac Computer Monitor
- > Stimuli™ Software Package
- > Remote Hand Controller
- > Pediatric Suite of Octotypes -  
    Tumbling E's and Tumbling C's
- > Cartoons Feature Package
- > Contrast Sensitivity Test



### Stimuli™ Advanced Version 4.0 Features:

- > Progressive, single line, single optotype and equal size charts, with random refresh, columns, isolation, pointer and contour bars.
- > Imperial or Metric Units for refraction calculation
- > Direct or Mirrored Installation
- > Multiple on-screen legends, including audible legends
- > Fixation Dot with options including reversed, animated and disappearing dot game.
- > Highest resolution Pixel Accurate™ optotypes
- > Pre-programmed randomized optotype
- > Sleek, all-in-one solution
- > Contrast Sensitivity test suite
- > Wireless remote optotype selector and pointer
- > Specialized Pediatric Test Suite

Feature Chart		Advanced
Sloan	<b>A B C D</b>	•
Snellen	<b>C D E F</b>	•
Numerals	<b>2 3 4 5</b>	•
Tumbling E	<b>W E M E</b>	•
Landolt C	<b>O C O O</b>	•
HOTV	<b>H O T V</b>	•
Ped Shapes		•
AC Shapes		•
T-Rings		•
Pointer		•
Letters		•
Duochrome		•
Fixation		•
Astigmatism Dial		•
Contrast		•
One Clip		•
Multi Clip		•
Animal Parade		•
Cartoons		•

### Exclusive Test Selection:

- > Astigmatism Dials
- > Contrast Sensitivity
- > Sizable Four Dots Test
- > Animal Parade
- > Selection of Optotypes



### Programmable Handset:



The Superior Practice.