


Veatch Manual

VISION TESTER



Focus,
To create delicacy

Each part manifests craftsmanship, bringing you splendid experience in remarkable ways.

Complete Assistant Glasses


Built-in multiple assistant glasses with comprehensive visual function examination.

First-class Feel

Precise socket design with clear location feedback of every movement.

Advanced Coating

Innovative lens coating, to reduce reflection and improve visual comfort degree.



Product Features

- Inner assistant glasses with cross optical lens and accurate color filter lens.
- Unique butterfly appearance, dazzling handle wheel, exquisite workmanship and comfortable feel.
- All optical lenses adopt advanced coating process.
- Accurate location when overturning cross cylinders.
- Comprehensive visual function examination and precise and convenient measurement.

Product Parameters

Sphere lens	-19.00m ⁻¹ ~+16.75m ⁻¹
	Step0.25m ⁻¹ (being 0.12m ⁻¹ when 0.12m ⁻¹ auxiliary lens is used)
Cylinder lens	0~-6.00m ⁻¹ (being 0~-8.00m ⁻¹ when additional lens is used)
	Step 0.25m ⁻¹ (being 0.12m ⁻¹ when additional lens is used)
Cylinder astigmatism axis angle	0°~180°, Step5°
Lens diameter	50~75mm, Step1mm
Rotating prism lens	0~20cm/m, Step1cm/m
Prism base	0°~360°
Cross cylinder	±0.25m ⁻¹
Lens pupil distance	13.75mm
Size and weight	335(W)×90(D)×320(H)mm/4.5kg