

4200

Pneumatic Stools

800.735.0357
haag-streit-usa.com

Backrest Adjusts Horizontally and Vertically on Model 4246

Model 4240
Seat Substrate 14-Ply Birch

High Density Polyfoam Seats

Lift Lever for Pneumatic Assist with Raising Seat

Exceptional Stability from Five-legged Design



4200

Pneumatic Stools



800.735.0357
haag-streit-usa.com

- ◆ **Finished for Durability**
 - Select from 14, easy-to-maintain vinyl colors. Other vinyl colors available by request.
 - All exposed metal, with exception of the aluminum base, is chrome-plated steel.
 - Base and legs finished in tough, black epoxy enamel.
 - Hooded dual-wheel casters.
- ◆ **Optional Foot Ring**

Order a foot ring with your stool, or at a later date.
- ◆ **Tech Specs**

Shipping weights: 110 (19 lbs.). 112 (26 lbs.). 4240 (23lbs.). 4246 (26 lbs.).

	4240	4246
Seat adjustment	Pneumatic Assist	Pneumatic Assist
Seat diameter	15"	15"
Seat construction	2" polyfoam	2" polyfoam
Seat swivels	360°	360°
Five legs for extra stability	Yes	Yes
Height range	19.50-27"	19.50-27"